

Points for Attention

1. To move this instrument, it is necessary to lift it up and put it down gently; in case of failure, please do not knock hard on it so as to prevent malfunction caused by the impact.
2. To move this instrument, it is necessary to keep it horizontal so as to avoid vibration and impact.
3. Do not squeeze the instrument body and do not place the heavy stuff onto the instrument body so as to prevent shell deformation, damage or malfunction.
4. It is necessary to select the corresponding working power supply in accordance with the marks on the instrument body.
5. It is not allowed to use this instrument in the bathroom.
6. To prevent the lightning accident, when it thunders, please turn off the power switch and pull out the power plug.
7. It is not allowed to inject any fluid into the instrument body so as to prevent the instrument from being damaged due to the short circuit.
8. If the instrument will not be used for a long time, it is necessary to pull out the power plug from the socket.

Technical parameters

1. Operating voltage AC 220 V--AC 240 V
2. Maximum power of the complete machine 25 W
3. Packing size 27*10*27(L*H*W) CM inner size

Version 20090808

USER MANUAL

Microcomputer Spot-removing Instrument



Dear customers;

We appreciate your order this product! Please go through this manual before use in order to ensure the maximum performance of the instrument and the safety in use. Please keep properly for further requires.


Functions

Microcomputer touch-control instant spot-removing instrument is a high-tech beauty instrument which generates the high-frequency current through the computer system, and removes the freckles, senile plaques, and polypi etc for an instant, and realizes an efficient and safe medical beauty result without bleeding and pains and any side effect.

Characteristics



1. It adopts the latest technology and full computer control and is easy and fast.
2. It utilizes adjustable sweep spot power output of 10 grades and is with smart and reliable treatment.
3. It has a unique sine wave pulse system and provides a variety of treatments.

Panel Operating Instructions



- | | |
|-------------------|---------------------------|
| 1. Time minus | 2. Energy minus |
| 3. Power on/off | 4. Time plus |
| 5. Energy minus | 6. Standby indicator lamp |
| 7. Switch on lamp | 8. Time display |

Operating Method

1. Plug in the power socket and turn on the power supply.
2. Hear the startup sound “Toot” and time display is 30; energy is zero, the indicator lamp is green and the machine is in a standby mode.
3. Press the “ON / OFF” button, the red starting up indicator is on, the machine starts and the time second point flashes.
4. Press the  button that is of spot-removing minus, press the  button that is of spot-removing energy plus while the energy lamp is increased and decreased by the energy size of 10 grades.
5. There is no over-voltage when taking metal for a test.
6. During operation, the electric current shall be appropriate because too high a current will cause a scorch but too weak a current will affect the operation time and result. It is advisable to let the electrical needle to gently touch the spot and move back and forth during spot removing. It can be clearly observed that the melanin becomes scorch and falls off after the scanning of electrical needle. If the spot is small and shallow in area, it can be removed only at a time. If the spot is large and has a deep pigment, then it may need to be treated twice to three times. During operation of electrical needle, the beautician shall pay attention to keep balanced and not wobble hands, so as to avoid burns and uneven treatment, without affecting the recovery of wound and leading to scars.

Removal of nevus

During removal of nevus, it is necessary to differentiate the progress of nevus. Some nevi grow on the skin surface but some may be rooted in the subcutaneous tissue so it is not advisable to make the electrical needle lacerate the surface of the nevus at a time, you had best burn it off in the form of small square. If the nevus is relatively deep, upon completion of treatment, it is advisable avoid irritation and apply some nutrition ointment to the wounds 23 days before recovery, so as to avoid the formation of pits on the skin. Don't make the wound subject to irritation when applying the ointment.

Removal of granulation

Apply the electrical needle end to gently touch the root of granulation. Apply the laser electrical needle to burn the root of granulation to death, and the granulation will automatically fall off about after one week. Or apply the laser electrical needle to remove the granulation generated at a time. But if the granulation is too large, it is necessary to remove it for a number of times.

Notice for use

1. If the beautician wear the sterile gloves, she or he shall first clean and sterilize her or his hands (by means of alcohol).
2. Fix the sterilized needle to the handle.
3. Clean the treated part of the customer with alcohol.
4. For the operation part of customer, wipe the part to be treated about 3 minutes with anaesthetic, and then start the operation of spot removing.