

Charging Method

1 Charging line connection USB interface

The battery is not fully filled when the product leave the factory. Please charge the battery fully before use

2 The USB interface of Computer, adapter (sold separately) and mobile phone could be used to recharge

LED lights flicker continuously in the process of charging, and will stop flicker when fully filled. Please confirm the state of the lithium battery before charging. Please tear the protective film down before the first use.



Notice

Please charge the device, if it unable to boot after long press the power button.

Please charge the device, if the device makes a sound "bi-bi-bi", the battery power shortage and the device will turn off automatically.

The product will be in standby state automatically after use ten minutes and will turn off after one minute if not used.

Fully charged needs about 150 minutes (First use, need to charge more than five hours.)

The device can be used about 80 minutes after fully charged.

Please charge the device in time before power off.

Instructions

Cleansing

- 1 Press the power button for 2 seconds to turn on device
- 2 Press the button one time to select the Cleansing mode
Cleansing mode, with color light function at the same time. Holding both sides of the conductor to enter the +ion mode. It could be better if clean the face with cleanser before start mode.
- 3 Smear the derived liquid on the face
- 4 As shown in Figure 2, clean face with the sonic head.



Moisturizing

- 1 Press the power button for 2 seconds to turn on device
- 2 Press the button two time to select the Moisturizing mode
Moisturizing mode, with color light function at the same time. Holding both sides of the conductor to enter the -ion mode. It could be better if use mask before start this mode.
- 3 Smear the beauty liquid on the face
- 4 As shown in Figure 2, leading-in with the sonic head



Usage Method

Lifting

- 1 Press the power button for 2 seconds to turn on device
- 2 Press the button three time to select the Moisturizing mode
Lifting mode, with color light function at the same time. It could be better to use with essence together.
- 3 Smear essence on the face
- 4 As shown in Figure 2, Lifting with the sonic head



Cleaning & Maintenance

- 1 Please make sure that the device is in power off state
- 2 Please clean up the device with cotton pads or soft fabric
DO NOT use diluent, solvent and other corrosive industrial products to clean sonic head in order to avoid damage
- 3 Please put the device in a dry and cool place
DO NOT put the device in the wet place, such as the bathroom
DO NOT put the device in a place where children can touch it

Sonic Skin Scrubber



Use Manual

Charging Method



Other Accessories ● Products Manual ● USB Charging Line

Product Introduction

Ultrasound is a high frequency mechanical vibration wave, which is widely used in medical and cosmetic fields. The mechanical wave has mechanical, warm and catalysis function. It can produce gentle resonance effect on skin at low energy and high frequency vibration, enhance cell viability and permeability of all cortical tissues, dredge sweat glands and shrink pores.

Ultrasonic cleaning and beauty instrument using nanotechnology, molecular 24000 times per second jolt of intense friction, the vibration deep into the skin to intense friction between molecules, deep temperature, cell metabolism, accelerate microcirculation, promote metabolism of back flow. Through the transformation of the ion flow, increased permeability of the cell membrane, the membrane to produce intracellular protoplasmic micro flow, fatty acid and H2O fatty CO2 alienation into delicate volume decreases, at the same time with the targeted products, cosmetic effect is more significant.

Deep massage through high-frequency vibration can improve microcirculation during fat movement, and make fat cells burn and lyse to metabolize products, so as to achieve clean effect. Using a special shovel type probe, can effectively penetrate into skin under 4-7 cm, let beauty liquid become thin, through supersonic vibration atomization, the dirt and excess sebum to emulsification, jam cosmetics residue, in the pores of sweat and emulsion, with a splash eliminated, the cleaning effect more quickly and effectively. This is the same as a supersonic wash that cleanse precision machinery, and this kind of phenomenon will occur on the skin. To remove the cuticle, fat particles, acne and so on, so that the nutrients of your maintenance products will no longer waste, so that the skin can be better absorbed.

Product Features

- 1 Deep cleaning and scavenging aged horniness
- 2 Pigment decomposition and desalination spot
- 3 Deep rhytidectomy and enhance skin elasticity
- 4 Promoting blood circulation and lymphatic metabolism
- 5 Enhance the regenerative ability and permeability of skin cells
- 6 Skin massage, nutrition introduction and variegated derivation
- 7 No damage and side effects on skin
- 8 Elegant appearance, easy to use, operate and carry

Product Parameters

Product Name: Sonic Skin Scrubber
Power Supply: AC100V-240V 50HZ/60HZ
Output Voltage: 5V/600mA
Charging Mode: Rechargeable (lithium battery)
Net Weight: About 99g
Power: 2W
Size: 173*55*15.5MM

Common Issues

| Issues | Cause analysis | Solution |
|----------------------------------|--|---|
| Power supply cannot be connected | <ul style="list-style-type: none"> Insufficient charge The battery power shortage | Please refer to the instruction "charging method" |
| Shutdown in use | <ul style="list-style-type: none"> Automatic shutdown function (auto-shut-off function after about 20 minutes) The battery power shortage (the device makes a sound "bi-bi-bi") | Normal function, non fault Please refer to the instruction "Charging Method" |
| Unable to recharge | <ul style="list-style-type: none"> USB interface is not connected with rechargeable device Lithium battery charging device for mobile phone used in charging The charging device failure The adapter (sold separately) fault | Please connect the USB charging line Please replace the charging device Please restart the charging device Please use other charging devices |