

SURGICAL FACE MASK – 3 PLY AEROFIT WHITE

ART NR 10001 (S-SERIES)

Product Data Sheet

Product Details

- One Size
- Earloop, Type IIR.
- Colour: White
- Hypoallergenic
- Non-sterile
- Single use / non-washable / non-reusable



Material Compositions

- Nose bar –Single core nose bar - adjustable
- Outer layer –Spunbond Non-wovens
- Middle layer – Meltblown filter layer
- Inner layer – Spunbond soft fiber layer.
- Earloop – flat elastic earloop -adjustable
- Material:
 - Odor free
 - No animal tissues
 - No human blood derivatives
 - No natural rubber latex
 - No medicinal substances or harmful chemicals (Health Care without Harm)
 - No Polyvinyl Chloride
 - No nickel

Country of Origin

- Made in Malaysia.

Shelf life and Storage

- 3 years from the date of manufacture.
- Use By Date on the box in MM/YYYY format.

Intended Usage

- Healthcare settings during infectious disease outbreaks.
- Industrial Settings & construction sites.
- Urban environments with high pollution levels.
- Always follow User Instructions and use in manner as indicated.

Packing Level	Quantity	GSI code	WxLxH(mm)	Art nr
Dispenser Box	50	9551019500049	175 x 95	10001

Testing Protocols and Results

Characteristics	Testing Protocol	Requirement	Product performance
Bacterial filtration efficiency (BFE)	BS EN 14683:2019	≥98%	99.9%
Differential pressure	BS EN 14683:2019	<6 mm H ₂ O/cm ²	1.726 mm H ₂ O/cm ²
Microbial cleanliness	BS EN 14683:2019	≤30 cfu/g	9 cfu/g
Particle filtration efficiency (PFE)	BS EN 14683:2019	≥98%	98.5%
Splash resistance pressure	BS EN 14683:2019	≥16 (kPa)	≥16 (kPa)

Applicable Regulation

EU directive MDR 2017/745 - class I (non-sterile), EN 14683 type II R

Dimensions

Dimension Text	Dimension Value (mm)
Length	95
Width unfolded	150
Width folded	175
Length ear-loop (unstretched)	170
Nose clip	105


Storage Instruction

- Product should be stored under normal storage conditions.

Front of mask – identifiable via Empro logo



Inside of mask – identifiable by reverse Empro logo



For More Information

Technical Assistance: (+60) 13- 3191919

Email: Marketing@empro.my

Hours of Operation: M-F 9 am – 6 pm

EMPRO Protective products are occupational use or private use only.