

# SURGICAL FACE MASK – 3 PLY AEROFIT GREY

ART NR 10005 (G-SERIES)

Product Data Sheet

## Product Details

- One Size
- Earloop, Type IIR.
- Colour: Grey
- Hypoallergenic
- Non-sterile
- Single use / non-washable / non-reusable

## Material Compositions

- Nose bar –Single core nose bar - adjustable
- Outer layer –Spunbond Non-wovens
- Middle layer – Meltblown filter layer
- Inner layer – Spunbond soft fiber layer.
- Earloop – flat elastic earloop -adjustable
- Material:
  - Odor free
  - No animal tissues
  - No human blood derivatives
  - No natural rubber latex
  - No medicinal substances or harmful chemicals (Health Care without Harm)
  - No Polyvinyl Chloride
  - No nickel



## Country of Origin

- Made in Malaysia.

## Shelf life and Storage

- 3 years from the date of manufacture.
- Use By Date on the box in MM/YYYY format.

## Intended Usage

- Healthcare settings during infectious disease outbreaks.
- Industrial Settings & construction sites.
- Urban environments with high pollution levels.
- Always follow User Instructions and use in manner as indicated.

Packing Level	Quantity	GSI code	WxLxH(mm)	Art nr
Dispenser Box	50	9551019500155	175 x 95	10001

## Testing Protocols and Results

Characteristics	Testing Protocol	Requirement	Product performance
Bacterial filtration efficiency (BFE)	BS EN 14683:2019	≥98%	99.9%
Differential pressure	BS EN 14683:2019	<6 mm H <sub>2</sub> O/cm <sup>2</sup>	1.726 mm H <sub>2</sub> O/cm <sup>2</sup>
Microbial cleanliness	BS EN 14683:2019	≤30 cfu/g	9 cfu/g
Particle filtration efficiency (PFE)	BS EN 14683:2019	≥98%	98.5%
Splash resistance pressure	BS EN 14683:2019	≥16 (kPa)	≥16 (kPa)

## Applicable Regulation

EU directive MDR 2017/745 - class I (non-sterile), EN 14683 type II R

## Dimensions

Dimension Text	Dimension Value (mm)
Length	<b>95</b>
Width unfolded	<b>150</b>
Width folded	<b>175</b>
Length ear-loop (unstretched)	<b>170</b>
Nose clip	<b>105</b>

## Storage Instruction

- Product should be stored under normal storage conditions.

**For More Information**

**Technical Assistance: (+60) 13- 3191919**

**Email: [Marketing@empro.my](mailto:Marketing@empro.my)**

**Hours of Operation: M-F 9 am – 6 pm**

**EMPRO Protective products are occupational use or private use only.**