

ANAGEL® Ultrasound Gel

Ultrasound gel for medical diagnosis and therapy.

Features

- Acoustically correct for the broad range of frequencies used
- Clear Gel
- Odourless
- Completely aqueous, will not stain clothing or damage transducers
- Bacteriostatic
- Non-sensitising non-irritating
- No formaldehyde or salt

Standard sizes: 20ml, 60ml, 250ml, 1000ml and 5L

Expiry date: 3 years from date of manufacture



SKU	Size	Case Size
UGEL20	20ml	800
UGEL60	60ml	144
UGEL250	250ml	40
UGEL1000	1000ml	12
UGEL5000	5000ml	4

Technical Specification:

ANALYSIS PARAMETERS:

RESULTS

REFERENCE METHOD

Physical and Chemical Analysis

Colour	Colourless	-
Smell	Characteristic odour	-
-		
Appearance	Gel	-
Density (20 C)	0,98364 gr/cm ³	ASTM D 891-89
Viscosity (25 C)	1,180,000 Cp	Brookfield DV-
II+Viscometer		
Ph(20 C)	6,75	Inolab Ph Meter