Technical Data Sheet

Anagel Arnica Gel v1.0

ARGEL250 – 250ml Bottle ARGEL500 – 500ml Bottle



Product:	Date:
Anagel Arnica Gel	February 2016
References:	Specification N.:
ARGEL100 - 100ml Tube	ETP/144/8

TECHNICAL SPECIFICATIONS

Composition (INCI)

Aqua (Water), Alcohol Denat, Glycerin, PEG-8, Polysorbate-20, Propylene Glycol, Carbomer, Triethanolamine, Lavandula Officinalis Flower Oil, Diazolidinyl Urea, Menthol, Arnica Montana Flower Extract, Camphor, CI 19140, CI 14720, Benzyl Alcohol, Limonene, Linalool, Coumarin.

Physical-chemical properties

- Colour, Odour, Appearance: Transparent yellow gel, characteristic odour to its raw materials (menthol, camphor).
- Viscosity (22°C) :100.000cps 160.000cps
- pH (22°C): 6,6 7,6
- Density (22°C): 0,95 1,00g/ml

Microbiological properties

- Total aerobes: <100ufc/ml
- Moulds and yeasts: <10ufc/ml
- Pathogens (Staphylococcus aureus, Candida albicans, Pseudomonas aeruginosa, Escherichia coli): Absence

PAO and expiry date

- Period after opening: 12 months
- Expiration: 3 years upon the manufacturing date.

Packaging

- ARGEL100: 100ml transparent PE tube with dosing cap
- ARGEL250: 250ml transparent PET bottle with dosing pump
- ARGEL500: 500ml transparent PET bottle with dosing pump

PRODUCT CHARACTERISTICS

- Gel with Arnica extract (Arnica Montana), Lavender essential oil (Lavandula Officinalis), Menthol and Camphor.
- Thanks to the beneficial properties of Arnica, the product provides wellness to your skin.
- Relieves, revitalizes and soothes.



Technical Data Sheet

Anagel Arnica Gel v1.0



INSTRUCTIONS FOR USE

Indications

Massage gently onto the desired area. This product can be applied in the entire body.

Warnings

For external use only.

Do not swallow.

Avoid contact with eyes, mucous membranes and wounds. Seek medical advice before using in children under 3 years of age, or if you are pregnant or breastfeeding.

REGULATORY INFORMATION

Ana Wiz Ltd guarantees that this product is in conformity with the European Regulation 1223/2009 for Cosmetic Products and that it has been manufactured following the directives of the Quality Assurance System based on ISO 22716.

CPNP Reference: 2109866