

Biological compatibility test

File No.: HL/CE-001-05

Version No.: B/6

Controlled state:

Issue number:

Prepared/date by: 周民 2017-12-26

Checked/date by: 俞素琼 2017-12-26

Approved/date by: 俞瑞兴 2017-12-26

File name: Biological compatibility test

File No.: HL/CE-001-05

The acupuncture needle manufacturing materials must have good biocompatibility, suitable for

Suzhou HuaLun Medical Appliance Co., Ltd.

version number: B/6

CE technical documents

effective date: 2017.12.26

medical, acupuncture needle or other material selection of SUS304 austenitic stainless steel wire, which is in strict accordance with the provisions of the national standard "GB2024-94" of the acupuncture needle and the choice of GB/T4240 "stainless steel wire" in SUS304 or other austenitic stainless steel wire to ensure the use of requirements, factory selection for acupuncture needle with SUS304 stainless steel wire, after years of use in full compliance with the requirements of the use of acupuncture needles, the material in accordance with the requirements of ISO7153-1:2000 with biological compatibility, no side effects on people.

As described above, the materials used in our company for acupuncture needles totally follow the requirements of china standards. They show good biocompatibility, ensure that the patient and the needle doctor do not cause severe toxicity or allergic reactions due to the manufacture of the material.

Enclosure:

Biological compatibility evaluation report

Stainless steel wire material report